

Fentanyl Resistance Statement

Eagle Protect, PBC hereby confirms that the product/s listed below have been independently tested by third-party laboratories for resistance to permeation by Fentanyl, as per ASTM D6978.

Products tested: Visible Nitrile Gloves – Powder Free, Black (6mil)
Blue 52 Nitrile Gloves – Powder Free, Blue (5mil)
Derma2 Nitrile Examination Gloves – Powder Free, Blue or Black (3mil)
Sensitive Nitrile Examination Gloves – Powder Free, Blue (4mil)

Results (applicable to all tested products):

Test conducted: Permeation testing per ASTM D6978

Test Drug and Concentration: Fentanyl Citrate Injection (100 mcg/2ml)

Test period: 240 minutes (4hours)

Test results: No breakthrough time up to 240 minutes.

Other observations: Slight swelling, no degradation.

This confirmation shall continue until revoked in writing by Eagle Protect, PBC. Conformance documentation is proprietary to Eagle Protect, PBC and its Suppliers, but may be provided for audit purposes on request, or redacted copies provided. Laboratory test reports have been shared by the product manufacturers with Eagle Protect, PBC, supplier of the product/s listed above.



Kath Row
Quality/Compliance
Eagle Protect, PBC



Eagle Protect PBC
800 384 3905
info@eagleprotect.com
eagleprotect.com

3079 Harrison Ave #21
South Lake Tahoe
California 96150
USA